

Quantiferon-TB Gold Plus: Results and Interpretation

Nil	TB1 – Nil	TB2 – Nil	Mitogen – Nil	Qualitative Result	Interpretation
≤8.0 IU/mL	≥0.35 IU/mL and >25% of Nil	Any	Any	Positive	<i>M. tuberculosis</i> infection likely
	Any	≥0.35 IU/mL and >25% of Nil			
	≤0.35 IU/mL OR ≥0.35 IU/mL and <25% of Nil		≥0.5 IU/mL	Negative	<i>M. tuberculosis</i> infection NOT likely
			<0.5 IU/mL	Indeterminate	Likelihood of <i>M. tuberculosis</i> infection cannot be determined
>8.0 IU/mL	Any				

QFT-Plus is considered positive if either one or both Tb Ag tubes are positive