



	<p><b>Copan eSwabs</b> Bacterial culture</p> <p><b>White: regular</b> Routine culture &amp; Rapid Strep</p> <p><b>Green: minitip</b> ear/eye/mid-meatus/SCAN/etc.</p> <p><b>Blue: flexible minitip</b> Pertussis</p>
	<p><b>UTM &amp; FLOQ swab</b> Viral PCR testing</p>
	<p><b>Adenovirus PCR Eye Swab Collection Kit</b></p>



## Sterile Petri dish

Dermatophyte culture – hair, nail, skin scrapings



## Sterile Cup

Body fluids, Tissues, Urine,  
Stool – C diff or Rotavirus



## Fecal Swab – culture AND Enterohemorrhagic *E. coli*

Feces



## Rectal Fecal Swab – culture ONLY

Does NOT include Enterohemorrhagic  
*E. coli* antigen test

Collar swab for Rectal sampling



## O&P Kit –

O&P Crypto/Giardia Antigen,  
O&P Foreign Travel

10% formalin & Zinc PVA



## Aptima Collection Kit

Trichomonas PCR  
GC & Chlamydia PCR\*

\*Not used for SAFE Exam/SCAN



## Chlamydia Collection Kit

GC & Chlamydia NAAT\*

\*Use this for SAFE Exam/SCAN



## Urogenital Wet Mount swab

Must add 1 mL saline