

	<p>Copan eSwabs Bacterial culture</p> <p>White: regular Routine culture & Rapid Strep</p> <p>Green: minitip ear/eye/mid-meatus/SCAN/etc.</p> <p>Blue: flexible minitip Pertussis</p>
	<p>UTM & FLOQ swab Viral PCR testing</p>
	<p>Adenovirus PCR Eye Swab Collection Kit</p>



Fecal Swab – culture AND Enterohemorrhagic *E. coli*

Feces



Rectal Fecal Swab – culture ONLY

Does NOT include Enterohemorrhagic *E. coli* antigen test

Collar swab for Rectal sampling



O&P Kit –

O&P Crypto/Giardia Antigen, O&P Foreign Travel

10% formalin & Zinc PVA



Chlamydia Collection Kit

GC & Chlamydia NAAT




Urogenital Wet Mount swab

Must add 1 mL saline



Sterile Petri dish

Dermatophyte culture – hair, nail, skin scrapings



Sterile Cup

Body fluids, Tissues, Urine, Stool –
C diff or Rotavirus