

## Copan eSwabs

Bacterial culture

**White: regular** Routine culture & Rapid Strep

**Green: minitip** ear/eye/mid-meatus/SCAN/etc.



## UTM & FLOQ swab

Viral PCR & Pertussis testing



## Adenovirus PCR Eye Swab Collection Kit









## Sterile Petri dish

Dermatophyte culture – hair, nail, skin scrapings



## Sterile Cup

Body fluids, Tissues, Urine, Stool – C diff or Rotavirus

	<p><b>Fecal Swab</b> – culture AND Enterohemorrhagic <i>E. coli</i></p> <p>Feces</p>
	<p><b>Rectal Fecal Swab</b> – culture ONLY</p> <p>Does NOT include Enterohemorrhagic <i>E. coli</i> antigen test</p> <p>Collar swab for Rectal sampling</p>
	<p><b>O&amp;P Kit</b> –</p> <p>O&amp;P Crypto/Giardia Antigen, O&amp;P Foreign Travel</p> <p>10% formalin &amp; Zinc PVA</p>
	<p><b>Aptima Collection Kit</b></p> <p>Trichomonas PCR GC &amp; Chlamydia PCR*</p> <p>*Not used for SAFE Exam/SCAN</p>
	<p><b>Chlamydia Collection Kit</b></p> <p>GC &amp; Chlamydia NAAT*</p> <p>*Use this for SAFE Exam/SCAN</p>
	<p><b>Urogenital Wet Mount swab</b></p> <p>Must add 1 mL saline</p>